

pro pen overview

Pro Pen offers both Microneedling and NanoInfusion options. 4 cartridges allow customization options to match the skin needs of each client. They can be paired with professional serums.



4 cartridge options for personalization

Grip for precision

Adjustment ring

Battery or cable operated to adjust to every environment. 2 batteries with each 2-hour charge = 4-hour charge

6 speeds: from 7500 to 8600 RPM

SMART technology remembers last speed setting





dermalogica PRO

cartridges

- Transparent color for more visibility
- Lock to prevent movement once plugged on pen
- Silicone spring prevents fluid backflow for maximum hygiene
- High grade sterilization method
- Stainless steel needles and plates



4 cartridge options allow for customization for every client.

				
	9-Pin	24-Pin	Nano Metal	Nano Silicon
	9 stainless-steel needles on square plate (5x5mm). This provides the most advanced form of exfoliation, smoothing and firming benefits.	24 stainless-steel needles tightly stacked together on a round plate (7.83mm diameter). This allows for once again a slightly stronger exfoliation and smoother skin.	Fine cones on a stainless-steel plate (5x5mm metal). The larger plate helps cover larger areas. It is designed for a slightly stronger exfoliation and helps with skin luminosity.	Fine cones on small silicon plate (2.87x2.87mm) . Very gentle, it is designed to delicately exfoliate when skin is somewhat sensitive.
material	Stainless steel	Stainless steel	Stainless steel	Silicon
plate dimension	5*5mm	7.83mm (diameter)	5*5mm	2.87*2.87mm

dermalogica PRO